

ImmunoTurbi

RF

QUANTITATIVE TURBIDIMETRIC TEST FOR DETERMINATION OF RHEUMATOID FACTORS

For *In-Vitro* Diagnostic Use Only

Store at 2°C to 8°C

OVERVIEW

A B Garrod in 1858 named the disease rheumatoid arthritis replacing the old terms arthritis deformans and rheumatic gout.

Rheumatoid arthritis, or RA, is an autoimmune and inflammatory disease, which means that your immune system attacks healthy cells in your body by mistake, causing inflammation (painful swelling) in the affected parts of the body. RA mainly attacks the joints, usually many joints at once.

Arthritis" literally means joint inflammation. Joints are places where two bones meet, such as your elbow or knee. There are many different types of arthritis with different causes and treatments. In some types, other organs, such as your eyes, heart, or skin, can also be affected.

INTENDED USE

Quantitative determination of rheumatoid factors in human serum/plasma. This test is for in vitro diagnostic use only.

PRINCIPLE

The Assay is based on antibody to Anti-Antibody (Antigen) reaction. The serum/plasma is added in latex solution. If the RF is present in the sample, it is attached to latex particles conjugated with human Gamma globulins which results in agglutination. The amount of agglutination is directly proportional to change in the absorbance in the sample. This quantity can be calibrated and measured photometrically.

CONTENTS OF KIT

1. Latex Diluent (R1)
2. Latex Reagent (R2)
3. Calibrator (R3)
4. Pack inserts

OPTIONAL MATERIAL REQUIRED

1. Semi-automatic or Fully Automated Biochemistry Analyzer
2. PPEs (Disposable Gloves, Mask, Safety Goggles, Apron)
3. Biohazard Dust Bin.
4. Test Tubes
5. Micropipettes

PRECAUTIONS /KIT STORAGE AND STABILITY

6. Please read all the information in this package insert before performing the test.
7. Do not use after the expiration date.
8. Store in between temperature 2°C to 8°C.
9. Do not use if damaged or leaked.
10. Do not open until you are ready to start the test.
11. Keep out of the reach of children.

WARNINGS

1. Do not reuse the tested reagent.
2. Follow the instruction to get accurate results.
3. Use appropriate personal protective equipment.
4. Dispose the leftover and used reagents and samples hygienically in biohazard waste.
5. Treat samples and reagent reaction volume as potentially infectious. Avoid contact with skin.
6. For in vitro diagnostic use. Not to be taken internally.

7. Do not mix the specimen sample or interchange the different specimen.
8. Do not use the reagents of other lots in combination with the kit.
9. Discard the remaining reagent in the kit.
10. The manufacturer and distributor of this product shall not be liable for any losses, liability, claims, costs or damages whether director cons sequential rising out for related to an incorrect diagnosis.

SPECIMEN COLLECTION

1. Fresh serum or plasma shall be used.
2. Serum shall be separated from freshly collected blood in plain tube or clot activator tube.
3. Freshly collected EDTA blood can be used to separate plasma for testing.
4. The samples stored at 2-8°C for 2 to 3 days can also be used.
5. Do not use highly haemolized or lipemic samples.

REAGENT PREPARATION

Working reagent Preparation: Mix the reagent R1 & R2 gently before use.

Mix the reagent R1 & R2 with the proportion as follows

8 ml Diluent + 2 ml Latex reagent (800µl+200µl).

Prepare the required amount of reagent only.

RF calibrator: Ready to use value mention on vial in IU/ml.

ONE POINT CALIBRATION (LINEAR RANGE UP TO 120 IU/ML)

RF Calibrator is ready to use.

TEST PROCEDURE

	CAL	SAMPLE
Sample	-	10 µl
Standard	10 µl	-
Reagent	1000 µl	1000 µl

1. Bring the kit component to room temperature before testing.
2. Take 400 µl of Diluent Reagent in cuvette and add 100 µl of latex reagent in it.
3. Put 5.0 µl of sample/control/calibrator in to the cuvette and mix.
4. Immediately read the absorbance of the reaction.
5. Read the absorbance of the reaction again after exactly 120 seconds.
6. Calculate concentration and interpret the results as per formula.

Wavelength	650 (600-650) nm
Cuvette	1 cm light path
Reaction Temperature	37 °C
Measurement	Against Distilled water
Reaction	2-point kinetics
Reaction Direction	Increasing
Sample / Reagent Ratio	1 : 100
Linearity	160 IU/mL

INTERPRETATION OF RESULTS

1. Concentration = $\frac{(A2-A1) \text{ Sample}}{(A2-A1) \text{ Calibrator}}$ x Concentration of calibrator
2. Reference Values: - Normal: 0 to 20 IU/ml, Positive: Above 20 IU/ml.

LINEARITY

The method is linear to a concentration of 6-160 IU/ml.

If the concentration exceeds this value, the sample should be diluted 1:5 with 0.9% saline solution and re-assayed.

LIMITATIONS

- 1. As with all diagnostic tests, the test result must always be correlated with clinical findings.
- 2. The results of test are to be interpreted within the epidemiological, clinical and therapeutic
- 3. context. When it seems indicated, reference correlation should be considered.
- 4. Any modification to the above procedure and / or uses of other reagents will invalidate
- 5. the test procedure.

QUALITY CONTROL

Positive and Negative controls are recommended to monitor the performance of test procedure, as well as a comparative pattern for a better results interpretation.

REFERENCE VALUE

Normal RF Value	0–20 IU/ml
Increased value	> 20 IU/ml












Each laboratory should establish its own reference range.

DISCLAIMER

The all precautions shall be taken to ensure the diagnostic ability and accuracy of this product. This product is utilized outside the control of manufacturer and distributors. The various factors including storage temperature, environmental conditions and procedure error may affect the results.

REFERENCES

- 1. Robert W Dorner et al. Clinica Chemica Acta 1987;167;1-21.
- 2. Vladimir Mule at al Scand J Rheumatology 1972;1;181-187.
- 3. Frederick Wolfe et al. Arthritis and Rheumatism 1991; 34:951-960

	In Vitro Diagnostic Use
	Manufacturer
	Manufacturing Date
	Expiry Date
	Lot Number
	Store at 2°C to 8°C
	Single Use
	Number of tests in the pack
	Do not use if pouch or kit damaged
	This side Up
	Read package insert before use